

## BIO PEN T6 – Wellness and Rehabilitation combined

### Introduction

Bio Pen T6, developed by Dr. Pen, is a professional grade device for skin tightening and lifting. It combines RF, EMS electroporation, and Blue light therapy to enhance skin texture, detoxify, minimize the appearance of enlarged pores, and reduce sun spots and acne scars. Using this advanced technology, professionals can assist their clients in achieving more youthful looking skin.



T6 is a beauty instrument that integrates four functions of radio frequency (RF), micro current (EMS), electric skin lifting and blue light therapy. It is suitable for beauty salons and home use. It has a powerful built-in intelligent control system. The adjustment of the three major functions of radio frequency, micro current and massage can realize customized classification of skin care for the forehead, eye area, face, jawline, neck, back, abdomen, legs, arms and other parts.

At the same time, it has an intelligent time control system, which can customize the working time to meet different needs; among them, the three modes of electric skin lifting can be superimposed with radio frequency, micro current and blue light to enhance the skin beautification effect.

### Advantages

- **Triple Photoelectric Power**

Featuring advanced skincare technology: Radio Frequency(RF), Electrical Muscle Stimulation (EMS), and Blue Light therapy to help you combat aging skin issues.

## Three lifting massage directions

Partition skin care & Massage lifting



### Clockwise lifting

It is often used for skin care of the obvious wrinkles and nasolabial folds on the right face, which is opposite to the growth direction of the pores to achieve an upward lifting trend.



### Anticlockwise lifting

It is often used for skin care of the obvious wrinkles and nasolabial folds on the left face, which is opposite to the growth direction of the pores to achieve an upward lifting trend.



### Both clockwise & anticlockwise soothing

It is often used in the steps of skin care massage to improve the absorption of skin care products, relieve fatigue, lift and shape the skin.




## Electric Skin Lifting

### Blue Light

#### EMS

EMS boosts cell activity, promote muscle movement and restore skin elasticity. It can increase blood circulation of capillaries and promote the absorption of nutrients by skin tissue.




## Electric Skin Lifting

### Radio Frequency

#### Blue Light


Skin heating therapy can stimulate the regeneration of collagen in the dermis. It can be combined with electric skin lifting to perform large-scale skin care. At the same time, it can control the skin temperature to improve skin care efficiency, reduce wrinkles and tighten the skin.



## Electric Skin Lifting

### Blue Light

Blue light can reduce inflammation and pigmentation, enhance the beauty effect, combined with Radio Frequency, EMS and Rotation, it can achieve twice the result with half the effort.



- **Customized Treatment**

Easily personalize each treatment session to suit your skin needs with the help of Customized Operation Assistance System and Built-in Time Control System.

- **Safe and Certified**

The Bio Pen T6 is safe, effective, and non-invasive proven with three certification for your peace of mind: CE-LVD, RoHS, CE-EMC.





## Custom Settings

In addition to the default values, RF intensity, EMS intensity, Rotation speed and working time can be reset according to each individual's needs to customize the process and obtain the ideal skin care effect.

